

Diabetes Technology Expert Program

CareSens Air

1. Introduction

Welcome to this module on the CareSens Air CGM. In this series of videos, we'll examine the ins and outs of one of the newest CGMs on the market. Let's begin.

Measuring approximately 3.5 by 2 cm and weighing 4.5 grams, the CareSens Air CGM is placed on the back of your upper arm using an applicator and secured with sensor tape. After a 30-minute warm-up, the sensor monitors glucose levels in your interstitial fluid and delivers real-time readings to the CareSens Air smartphone app or the CareSens Air receiver via Bluetooth every five minutes. The sensor lasts for up to 15 days. It is intended for people with diabetes 18 years or older.

In our “Master Module on Glucose Sensors,” we covered CGM fundamentals, including when to trust sensor readings and how to manage issues such as disconnections and skin reactions. This module centers on the start-up education for the CareSens Air CGM and adds advanced tips and tricks.

This module will help both healthcare providers and individuals with diabetes get the most from the CareSens Air system. But before we start, please note that these videos do not serve as medical advice. Always consult with your healthcare provider before modifying your therapy in any way.

These videos were produced by Dr. Inge Van Boxelaer, the founder of Diabetotech and an endocrinologist at the St-Lucas Hospital in Ghent, Belgium. The content was reviewed by Maike Henssen, a diabetes educator from Bosman in the Netherlands. With their expertise and perspectives, you are guaranteed a balanced assessment of the CareSens Air CGM.

Be assured, all the information provided here is completely independent from i-SENS, the sensor manufacturer, ensuring an unbiased evaluation of the strengths and limitations of the CareSens Air.

By the end of this module, you will have a thorough grasp of the CareSens Air CGM. We will discuss sensor placement and start-up, alarm settings, data sharing, device integrations, and best practices. As a healthcare provider, you will be equipped to confidently support users of the CareSens Air, while users themselves will learn practical strategies to optimize their daily experience with the device. Together, let's start this empowering journey. Best of luck!

2. Placing the sensor

Welcome to the first lesson of the CareSens Air CGM module. In this video, we'll walk you through placing and removing the CareSens Air sensor, along with key dos and don'ts. Let's dive in.

2.1 Supplies

To start, you'll need a CareSens Air kit, which consists of:

- An applicator with a built-in sensor;
- Sensor tape; and
- An alcohol swab, which may or may not be included in the box depending on regulatory approvals in your country.

2.2 Placing the sensor

Download the CareSens Air app, and create an i-SENS account and sign in. Scan the barcode on the sensor package label to enter the sensor information.

To prevent infection, wash your hands before attachment. Then open the CareSens Air sensor package. Take the applicator out of the package and set it on a flat, stable surface.

Choose the area where it will be attached. The sensor should be attached to the back of the upper arm, and the area should have enough fat to avoid muscles and bones. Clean the area thoroughly with an alcohol swab and let it dry completely.

Hold the safety cap down with one hand, and grab the applicator with the other hand to pull upward to remove the safety cap.

Place the applicator on the area. Be careful not to press the applicator onto the skin too hard or tilt it.

Press the blue release button on the applicator. Then remove the applicator. Make sure that the sensor has been properly attached. Gently press the white tape around the sensor to secure the placement.

2.3 Removing the sensor

The CareSens Air sensor can be used for up to 15 days. The sensor will be automatically disconnected when it expires.

If you wish to replace the sensor before it expires, you need to manually disconnect it. Tap the Settings icon at the bottom of the app home screen to manually disconnect the sensor. Tap Disconnect on the right side of Connection Management. The pop-up window will appear, asking you to enter the four-digit code to disconnect. Enter the four-digit code.

Once the sensor is successfully disconnected, the status will be shown as not connected.

Remove the sensor slowly from the edge of the sensor tape.

2.4 Do's & don'ts

Here are essential guidelines for placement and removal of the CareSens Air sensor.

- CareSens Air has not been evaluated or approved for the following individuals:
 - Infants and children under 18 years of age
 - Pregnant and breastfeeding women, and
 - Individuals that are on dialysis or critically ill.
- Do not use the sensor if the package is damaged or opened, or if the sensor itself is damaged or defective, as this may lead to infection.
- Do not remove the safety cap of the applicator until you are ready to attach the sensor.
- After you separate the safety cap from the applicator, be careful not to face it toward anybody.

- If a sensor part breaks or disconnects in the process of attaching the sensor, you must check whether the part has remained under the skin. If you experience inflammation, redness, swelling, or pain due to an infection at the site where the sensor was attached, seek assistance from your healthcare provider.

You've now learned the steps for placing and removing the CareSens Air sensor. Following these best practices helps ensure reliable glucose readings day after day, supporting confident treatment decisions.

3. Starting the sensor

Welcome to this lesson on starting the CareSens Air sensor. In this video, you'll learn how to start the sensor and navigate the app. Plus, we'll tackle practical tips that will allow you to make the most of your CareSens Air system. Let's jump right in.

3.1 Starting the sensor

To connect the sensor with the app, press the power button until it clicks. After the sensor is successfully connected, you can configure your alert settings. We'll dive deeper into alerts in the next lesson.

Sensor warm-up begins automatically and takes 30 minutes. Your glucose level will appear on the home screen once the sensor is ready.

3.2 The CareSens Air app

In the CareSens Air app, instant glucose data and glucose trend arrows are displayed on the home screen. Your dynamic glucose data will appear in different colors depending on where it is in your target glucose range. It will appear green if it is within, orange if it is high, and red if it is low.

Your glucose reading will be tagged as HIGH if it is greater than 500 mg/dL (27.8 mmol/L), and LOW if it is less than 40 mg/dL (2.2 mmol/L). If no glucose data has been received within the last 25 minutes, a dash sign will appear.

The Glucose Trend Arrow indicates the direction and the rate of change in blood glucose compared to the previous reading. It is displayed as an arrow on the right side of the latest glucose reading on the home screen.

There are seven trend arrows indicating: stable, slowly increasing, increasing, rapidly increasing, slowly decreasing, decreasing, and rapidly decreasing glucose changes. If there is insufficient data to determine direction and rate, the change is considered “Unknown,” and three dots (...) are displayed.

Tap the Glucose Statistics icon at the top left corner of the home screen to view the statistics of glucose levels in the past 24 hours. Tap the Screen Mode icon on the home screen to see your glucose trend in more detail. In the horizontal mode, you can see recorded events along with your glucose trends. Each event is displayed with an icon above the timeline. Tap the event icon for more details. You can see the 24-hour glucose trends and events. Zoom in with your fingers to see 6- or 12-hour trends.

If you need to reinstall the app due to changing your smart device or other reasons, you can download the previous data uploaded in the cloud server. After reinstalling the CareSens Air app and logging in with the previous account, Data Restoration page will pop up. To restore data while using the app, tap on the arrow on the top right corner of Previous Data.

You can add a variety of events related to glucose care in the Logbook page. Tap Logbook at the bottom of the home screen. Recorded event logs are displayed from the most recent events. To record a new event, tap the plus icon at the bottom right corner of the screen. The Event Record page will appear. Choose an event you want to record from Blood Glucose, Exercise, Insulin, Ketone, Meal, Medication, or Note. Fill in the details including the event time and date. Then tap Save. The event is recorded. You can find the recorded events from Logbook or Glucose Trends in horizontal mode.

You may choose to calibrate your glucose readings if you feel the need or think it does not reflect your physical condition, particularly when the difference between the sensor reading and a finger-prick test exceeds 20%. Remember to only calibrate your sensor if your glucose is stable and in range. CareSens Air supports a widget so you can see key info on your smart device's main screen. Use the widget to check your sensor connection status, glucose readings, and trend arrows, or to enter a calibration value. Android offers three widget layouts; iOS offers one.

You can enter a calibration value when the Calibration icon is blue. Tap the icon to open the Calibration screen. On the other hand, calibration is not available when the icon is gray.

Meanwhile, tap the Plus icon to open the "Add New Event" screen.

3.3 The CareSens Air receiver

If you cannot use the CareSens Air smartphone app or choose not to, you can use the CareSens Air Receiver instead. The receiver is operated using four buttons: a round Select/Next button, a Previous button, and Move Menu Up and Down buttons. Refer to the table for the other functions these buttons perform.

To turn on the CareSens Air Receiver, press the round button. Review the device terms, then tap Agree. Next, complete the initial setup by selecting the date format, entering the date, picking the time format, setting the time, and choosing the unit.

After reading the product description and safety information, connect the receiver to the sensor by entering the six-digit PIN found on the package label. If pairing is successful, you'll see a confirmation message. If the connection is cancelled, a cancellation pop-up will appear instead.

Once pairing is complete, you will be prompted to set threshold values for Very Low, Low, and High glucose alerts, as well as the Rapidly Changing Alert settings. You will also select your preferred alert type.

Once the 30-minute warm-up is complete, the receiver will take you to the Home Screen where you'll see the latest CGM value, the trend arrow, the unit, the CGM trend graph, the average blood glucose for the graph period shown, and the time remaining until the sensor expires.

To open the Menu from the Home Screen, press the round button. Menu options include Start New Sensor, Calibration, Events, Report, Information, Settings, and Shutdown. Selecting Events opens a sub-menu where you can enter event records for blood glucose, ketones, insulin, meals, oral medication, and exercise. You can also access the event log from the same section.

3.4 Do's & don'ts

Now, let's highlight some do's and don'ts for using the CareSens Air system.

- To connect the sensor with the app, ensure Bluetooth is enabled on your smartphone.
- For iOS users:
 - Enter your PIN code and tap "Pair" to finish connecting the sensor with the app.
 - Do not close the app after connecting the sensor. Otherwise, they will be disconnected.
 - You will not receive glucose data in real time if the app is closed in the background. Make sure to reopen the app once it is closed.
- If the sensor fails to warm up, contact the i-SENS customer service center or refer to the company website.
- A sensor in use with a logged-in account can be reconnected to the same device as long as the sensor remains valid. Once the sensor is connected to the app, it cannot be connected to another app. When reconnecting or changing devices, the sensor can be connected only to the app that it was first connected to.
- The glucose value shown in the iOS widget may differ from the most recent reading. Open the app to check the latest value.
- Keep the sensor and your smartphone within 6 meters or 20 feet of each other, without any obstacles such as walls or metal objects in between. The distance from the sensor to the smart device must be closer if there is any solid object between them. Otherwise, the connection may fail.
- The sensor is waterproof for up to 24 hours at a depth of 1 meter. You can wear it during a shower or bath or while swimming. However, it is recommended to attach waterproof tape to ensure safe use of the sensor when exposed to water for extended periods.
- If a skin care product such as sunscreen or insect repellent gets on the sensor, wipe it immediately with a clean cloth. Do not wash the sensor. Using an unsuitable solution could damage the device.
- Remove the sensor before an X-ray, MRI, or CT scan.
- Comply with your country's regulations for handling and managing medical waste when you dispose of the sensor.

We've now covered the CareSens Air system set-up, from sensor warm-up, app navigation and glucose calibration, to important reminders for proper use. As you support your patients in their first days with the CareSens Air, acknowledge that technology may sometimes come with operational bumps, such as sensor errors or connection issues. Should they arise, stay patient and persistent. Fully explore the available troubleshooting steps, and contact CareSens customer support as needed. With your knowledge and determination, you are well positioned to help your patients realize the full benefits of their CGM.

4. Setting alerts

Welcome to this lesson on setting alerts with the CareSens Air CGM. In this video, you'll learn about the types of alerts you can receive from the CareSens Air app and how to customize them. Let's get started.

4.1 Available alerts

The CareSens Air app sends alerts for:

- Glucose level (Very Low, Low, High)
- Rapidly changing glucose level
- Signal anomaly detection
- Signal loss
- Sensor replacement
- Sensor error
- Care provider connection

With the app running, the seven alert types appear as popup alerts. If the app is closed or your smart device screen is locked, you'll see the alerts as banners.

Glucose Level alerts occur when your readings fall below your Low or Very Low thresholds or above your High threshold. These alerts require treatment decisions.

You'll get a Rapidly Changing Glucose Level alert if your glucose changes faster than the rate you set. You may have to make treatment decisions in response to these alerts.

Signal Anomaly Detection alerts result from factors like improper sensor attachment or pressure on the sensor, which cause temporary inaccuracies in blood glucose readings. If you receive this alert, check that the sensor is properly attached, ensure nothing is pressing on it, and wait for readings to stabilize. Use a glucose meter if treatment decisions are needed during a signal anomaly.

A Signal Loss alert arises when the connection between your sensor and smart device is lost for 25 minutes or longer.

You'll receive Sensor Replacement alerts 5 days, 3 days, 1 day, and 1 hour before the sensor expires. Recall that the sensor can be used for a maximum of 15 days.

If an error in the sensor occurs during warm-up, 'Sensor Warmup Failure' will appear. If an error occurs after warmup has been completed, a 'Sensor Error' alert will occur. When there is a Sensor Error, the app disconnects from the sensor and can no longer be used. Remove the sensor from your arm as instructed.

Last but not least, you'll get a Care Provider Connection alert when a new care provider is registered in the app. We'll explore data sharing with care providers in the next lesson.

With the CareSens Air Receiver, you'll receive the same glucose and trend alerts that are generated in the CareSens Air app. Here's how those glucose alerts appear on the receiver.

In addition, the Receiver provides system alerts for Sensor Warmup,, Signal Loss, Sensor Expiration, Sensor Error, Abnormal Signal and Low Battery.

4.2 Setting alerts on the CareSens Air app

You can configure your alert settings right after you've successfully connected the sensor to your app, and right before sensor warm-up begins.

Later on, you can customize alerts for out-of-range glucose, rapidly changing glucose levels, system alerts, and text-to-speech. Start by tapping the Settings icon on the home screen.

You can customize Glucose Level alerts by setting the VERY LOW, LOW, and HIGH threshold values; the alert type; and whether or not to activate them. Switch the alerts on and off based on your preferences.

For Rapidly Changing Glucose Level alerts, you can choose a rapid-change threshold of above 2 or 3 mg/dL/min or above 0.1 or 0.15 mmol/l/min.

Tap System Alerts to change settings for Signal Loss, Sensor Replacement, Sensor Error, and Signal Anomaly Detection.

When you enable Mute Mode, alerts are silenced. Once you turn it off, the alerts return to your previously selected alert settings.

With Text-to-Speech on, you'll hear an additional voice notification whenever any of the seven alerts is triggered.

Finally, you can find out more details on alerts by tapping the bell-shaped icon in the upper right corner of the home screen. You can see the history of glucose alerts, sensor alerts, and care provider's connection alerts that you've received while using the app.

4.3 Setting alerts on the CareSens Air Receiver

From the Main Menu, go to Settings, then Alerts. Here you can customise the Very Low, Low, High, and Rapidly Changing Glucose alerts. You can also adjust system alerts, enable Mute mode, and reset your alert settings. Use the Up or Down arrow to select the alert you want to configure, then press the round button.

In particular, you can change thresholds for Very Low, Low, High, and Rapidly Changing alerts. By default, the Receiver sends these alerts with sound cues, but the alert type can be modified as well. Save your selections by pressing the round button.

From the Alerts sub-menu, select System to view the system alerts. Use the Up and Down buttons to navigate, then press the Round button to confirm your selection. Keep in mind that for system alerts, only the alert type can be adjusted.

You can also silence notifications by selecting Mute from the Alerts sub-menu. Confirm that you want to enable Mute mode by pressing the Round button. Keep in mind that Very Low, Signal Loss, Sensor Error, Sensor Warmup, and Sensor Expiration alerts are not affected by Mute mode.

Finally, you can restore all alerts to their default sound settings by selecting Reset Alerts.

4.4 Do's & don'ts

Keep the following tips in mind when setting alerts with the CareSens Air.

- To receive alerts from the app, you must enable notifications in your smart device's settings. Go to Settings on your device. Tap Notifications. A list of apps that send notifications will appear. Find "CareSens Air" on the Notifications screen and tap it. In the app details, tap Notifications and turn on "Show notifications."
- When customizing your alert settings, you can choose from Sound, Vibrate, Sound + Vibrate, or Mute. It may be difficult to distinguish CareSens Air alerts from notifications sent by your device or other apps when the app's alert type is set to sound or vibration. Adjust your alert settings to more easily recognizable options to help you make appropriate treatment decisions.
- Very Low Glucose, Signal Loss, Sensor Error, Sensor Warm-up, and Sensor Expiration alerts are not included in Mute Mode. These alerts must be muted manually.
- A pop-up alert will always be displayed, even if the smart device is in Mute or "Do not disturb" mode.
- If multiple alerts are triggered at the same time, the alerts will be delivered in the following order of priority: Signal Anomaly Detection alert > Glucose Level alert > Rapidly Changing Glucose Level alert.

We've now examined how the CareSens Air alert system works. By leveraging these customizable alerts, you transform potential challenges in your diabetes management into actionable moments. Ultimately, you're empowered to stay on top of your well-being.

5. Sharing data

Welcome to this lesson on sharing data with the CareSens Air CGM. This video will explain how you can share your glucose data with your healthcare provider and loved ones via the Sens365 app and Sens365 Web. We'll also go over the details of the available reports.

5.1 Sharing data with your care provider via the Sens365 app

Care providers can install the Sens365 app on their smart devices to receive and monitor your glucose data. To do this, you must first generate a six-digit code on the CareSens Air app.

Tap Settings at the bottom of the home screen. Tap Care Providers. Tap Connect to link your account to your care provider's account. Tap Generate Sharing Code. Choose the data you want to share with your care provider, then tap "Generate Sharing Code." Then share the generated code shown on the screen with your care provider.

Now let's enter the shared code on the Sens365 app. After running the Sens365 app on the care provider's smart device, tap the person icon at the bottom of the home screen. Tap Add. Enter the shared code provided by a sensor user and tap Ok. You can also enter the sharing code via the Settings at the bottom of the home screen. Tap Connections. Tap the Add Person icon at the upper right corner of the screen. Enter the shared code provided by a sensor user. Tap Ok to complete the process. Once completed, a notification pop-up will appear on your CareSens Air app.

You can also enter the code provided by your care provider on your CareSens Air app to connect with them. Open the CareSens Air app and tap Care Providers in the Settings, then tap Connect. Tap Enter Sharing Code. Enter the code from your care provider. If the code is valid, a care-provider confirmation alert will pop up on the screen. Make sure you are adding the right care provider. When the process is completed, a completion notification will appear.

You can temporarily stop sharing glucose data with your care provider or change items you share with them. At the bottom of the home screen, tap Settings, then tap Care Providers. On the Care Providers screen, tap the care provider whose sharing options you want to change.

Tap on Share Data to activate or deactivate sharing data. Tap the check boxes to stop sharing data on glucose reading and events or alerts.

If you no longer want to share data with a care provider, you can turn off sharing. On the Care Providers screen, tap the Delete icon next to the provider you want to stop sharing data with. Tap Disconnect on the pop-up to complete the process.

5.2 Sharing data with your care provider via the Sens365 Web

All data saved by the CareSens Air app on your smart device is automatically stored on the cloud server. To refresh the upload of your app data, tap Settings at the bottom of the home screen, then tap Upload.

You can view and manage your uploaded glucose data using the cloud-based platform Sens365 Web. Open your web browser and go to <https://sens365.i-sens.com>. Click Get Started on the start screen.

The first time you visit the site, you'll be prompted to choose your country of residence and language. Click Submit to confirm your selections.

Click create an account on the login screen if you don't have a Sens365 account yet. Enter and verify your email address. Accept the required terms, create your password and fill out all the required information. A "Registration completed successfully" pop-up will appear.

If you already have a Sens365 account, simply enter your registered email and password on the screen to log in.

On your first login to Sens365 Web, the user mode selection screen will be displayed. If you are a CareSens Air user, select Sensor User. If you are a care provider, select Care Provider.

Sensor users must click OK to confirm their country of residence and glucose reading unit. They will then be taken to the home screen.

To connect to a new care provider, click Connections in the navigation sidebar on the left. Click the Connect with Care Provider button on the right. Enter the sharing code of the care provider and click OK. Verify the name of the sharing group. After reading the agreement, check the box and click OK. Once completed, you can find the new care provider added on the list.

To disconnect from a care provider, click the Disconnect button, then Click Disconnect on the pop-up window.

Now, let's discuss Care Provider mode. When you log in, the screen displays group names and a list of sensor users in each group. A summary of each user's information is also displayed. You can drag a column to reposition it.

When clicking on the Dashboard button for a specific user, the Add to Dashboard pop-up appears. Assign a care provider to the user, then click OK. The dashboard will appear on the right side of the screen. You can open the dashboard at any time by clicking on the Dashboard button in the top-right corner.

Groups can be created or deleted. To create a group, click the three dots next to the group name. In the Group Settings pop-up, customize the group name and generate the group's sharing code, then click Save.

You can also add a new sensor user or invite one to the group. If you are adding a sensor user, enter the sharing code you received from them. If you are inviting a sensor user, confirm the group name, then send your sharing code to the user. You can copy the code or send email invitations directly. Your care provider sharing code remains the same unless you change it; you can generate a new sharing code as needed.

Moreover, you can invite another care provider to manage the user groups with you. To do this, click the three-dot menu next to the group name and select Invite Care Provider. Enter the care provider's email address, then choose Admin or Coordinator Authority. Select the acknowledgment checkbox and click Send.

5.3 Reports

In Care Provider mode, click the sensor user's name to access the Reports screen. By default, you'll be taken to the AGP report. On the left, you'll find Overview. On the right, you can access Comparison, Patterns, Trends, Statistics, and Events.

In the row above, you can adjust the time period you want to view.

The Overview page displays Time in Range, Average Glucose, GMI, Glucose Variability, and Standard Deviation. Daily Average Glucose, Insulin, Daily Time in Ranges, Time CGM Active, and Device Information are also shown. A customization option is available, allowing you to adjust the layout by adding, deleting, or moving blocks.

AGP is a standardized glucose report that includes Time in Range, Time CGM Active, Glucose Metrics, Ambulatory Glucose Profile, and Daily Glucose Profiles. You can check the percentage of time spent in the target glucose range using the TIR graph and estimate A1c using GMI.

On the Comparison page, two graphs are displayed based on the selected time periods. This makes it easy to see changes in glucose trends and improvements in glucose management by comparing the two graphs.

On the Patterns page, low, target, and high glucose ranges are marked on the graph for the selected time period. This helps you identify patterns by showing when low or high glucose levels occur and which events may be associated with glucose changes.

The Trends page shows daily and weekly glucose trends. You can view day-to-day fluctuations or see all daily trends together in a weekly graph.

The Statistics page shows glucose readings by day of the week and hour of the day. You can review daily and hourly average glucose statistics within the selected time period.

The Events page shows all event logs recorded in the CareSens Air app in chronological order.

Further up the screen, next to the sensor user's name, click Target Range Settings to adjust the target blood glucose range for daytime and night-time periods and tailor management to the user's goals.

In both Sensor User and Care Provider modes, you can download, print, and send reports, as well as export data. To download a report, click the downward-arrow icon in the top-right corner; the report will be saved as a PDF. To print, click the printer icon. To send the report, click the envelope icon; it will be sent as a PDF to the email address you specify. Finally, to export your data in CSV format, click the Export Data icon.

Meanwhile, on the CareSens Air Receiver, the Statistics option in the Main Menu lets you view reports for the past 7, 14, 30, or 90 days. Each report shows your average glucose, standard deviation (SD), and coefficient of variation (CV) for the selected period, along with your Time in Range.

5.4 Tips

Review these reminders before sharing your CareSens Air data:

- When connecting to a care provider via the CareSens app, the time limit for entering the generated sharing code is 5 minutes. If the code expires, create a new one.
- When you change your sharing options for a care provider, that provider is notified of the changes.
- Sens365 is a secondary viewer for the CareSens Air CGM System and is not intended to provide treatment decisions or to be used as a substitute for professional healthcare advice. Sens365 does not guarantee the real-time reception of glucose data of the CareSens Air CGM user.
- You cannot share data directly from the CareSens Air Receiver. To access data stored on the device, visit the Sens365 website, sign in, and upload the data from the receiver to the cloud server.

In this video, we tackled various options for sharing and uploading your CareSens Air data. With these features, you can work more closely with your care team, helping you make timely, informed decisions for better health.

6. Connecting other devices

Welcome to this lesson on connecting your CareSens Air CGM to other devices. In this video, you'll discover how to integrate your CGM with blood glucose meters and explore options with smartwatches, smart pens, and insulin pumps. Let's start.

6.1 Blood glucose meters

CareSens blood glucose meters with Bluetooth functionality can be linked to the CareSens Air app to automatically send glucose measurements. Tap Settings at the bottom of the home screen. Tap Blood Glucose Meter. Tap the Blood Glucose Meter icon at the top right corner of the screen to connect with a new glucose meter. Tap the blood glucose meter you wish to pair within the list of connectable devices. Enter the PIN code shown in the glucose meter display. Be sure to follow the specific instructions for your BGM model. Once pairing is complete, data from your previous blood glucose meter can be imported. Choose data you wish to import and tap Save.

6.2 Smartwatches

You can connect the CareSens Air app to a Galaxy watch or an Apple watch, provided the minimum system requirements are met. Visit the official CareSens Air website (<https://caresensair.com/content/compatibility>) to find smartwatches that have passed i-SENS's compatibility tests.

To connect your Galaxy Watch to your phone app, launch the Play Store app on your smartwatch. Search for "CareSens Air," then tap Install.

For iOS users, Search for "CareSens Air" in the App Store and tap Install. Launch the Watch app on the iPhone where the CareSens Air app is installed.

If the watch app is already installed, close and restart it. The CareSens Air app will appear in the AVAILABLE APPS list, and the watch app will install automatically on your Apple watch. Installation may take a few minutes to complete.

On the watch app screen, you'll see the status of the Bluetooth connection between your smart device and sensor; the date and time when the blood glucose level was last received; the most recently measured glucose level; the current rate of change in your glucose level compared with the previous measurement; and a graph displaying changes in your glucose level while the sensor is in use.

You can choose from the following four types of widgets for Galaxy smartwatches:

- Glucose trends
- Glucose level
- Gauge + glucose level, and
- Reception date and time + glucose level

For Apple watches, five widget options include different combinations of gauge, glucose level, reception date and time, chart data, trend arrows, and changes.

When the app is closed or your smartwatch is locked, your watch can still display glucose and sensor alerts. Just ensure that alerts are set to be sent to both your smart device and smartwatch in your device settings.

Take note that the smartwatch communicates with the mobile phone, not the sensor. If the smartwatch is not connected to the phone and the phone is not

connected to the sensor, the smartwatch will not receive alerts or glucose readings.

Finally, you can connect other smartwatches to your CareSens CGM via xDrip+. In the CareSens Air app, tap Settings, tap Manage Data & Connections, and turn on the switch next to xDrip+. Then, in xDrip+, go to Settings and Hardware Data Source, and select CareSens Air.

6.3 *Smart insulin pens*

At this time, the CareSens Air system does not support direct connectivity with smart insulin pens, nor can it be integrated with the Glooko Health app for this purpose. However, these integrations may be considered for future versions of CareSens Air.

6.4 *Smart insulin pumps*

The CareSens Air sensor does not integrate with commercial AID systems, although it is compatible with open-source AID systems such as AndroidAPS via xDrip. iPhone-based open-source systems do not yet support CareSens Air.

Be aware that using open-source apps like xDrip+ carries inherent risks and is not supported by i-SENS or approved by regulatory authorities. We provide this information transparently, though we don't recommend engaging in this practice without understanding the risks fully.

In this video, you learned how to connect your CareSens Air CGM with compatible blood glucose meters and smart devices. By putting these integrations to work, you can create a seamless setup that empowers more informed choices and tailored diabetes management.

7. Conclusion

Congratulations on completing this module on the CareSens Air! Learning how this CGM works is a significant step toward managing diabetes with greater confidence and effectiveness.

Before we end, we recommend that you watch this independently produced user video on placing the CareSens Air sensor:

https://www.instagram.com/reel/DAJM_47MCEY/?hl=en

It offers a practical, real-world demonstration that complements the theoretical knowledge you've gained in this course and can help you feel more comfortable applying what you've learned.

And always remember, you're not doing this alone. Seek support from your healthcare team, fellow individuals with diabetes, and online communities. Continue using technology as a tool to stay ahead and make informed decisions that support your health.

Once more, congratulations on finishing this module. May your knowledge and dedication keep moving you forward toward your goals.