

Kaleido

Homescreen and buttons

BASAL
Once you're pumping, your Basal section will show the basal rate you are currently receiving. If you need to change or review your basal profile, just navigate to **Basal** and select it by pressing the Enter button.

RESOURCE LEVELS
Your Resource Levels section shows your current Insulin on Board (IoB) on your home screen. If you'd like to check all your Resource Levels, navigate to Resource Levels/IoB and press Enter. Once you've done this, the Resource Levels screen will display your current battery level, insulin cartridge volume, IoB and time until cartridge replacement



BOLUS
Your Bolus section shows you the details of the last bolus you delivered – that's the date, time and the amount of insulin you got. When you need to deliver your next bolus, navigate to **Bolus** using the buttons and press Enter.

MENU
If you'd like to take a look at your History, upload data or check or update your personal information or system settings, you'll need to navigate to **Menu** and press Enter.

STATUS
To pause, stop and restart your pump, navigate to **Status**, press Enter and go from there.

BUTTONS
You can navigate around your Kaleido handset simply by using the up, down, left, right and Enter buttons.

Menu Map

